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Innovative Circular Dichroism and Circularly Polarized Luminescence Measurement Methods
(¹JASCO Co.) ○Yoshiro Kondo,¹

Circularly polarized luminescence (CPL) spectroscopy measures the difference in fluorescence intensity between left- and right-circularly polarized light. CPL spectroscopy has been used to analyze the structure of chiral molecules in their excited state. Recently, it has found an increasing range of applications in the analysis of molecules that emit circularly polarized light and can be used in 3D displays, and the number of articles focusing on CPL spectroscopy has increased dramatically.

Until now, the principal targets of CPL measurements have been solution samples. However, for practical device applications, it is also necessary to measure CPL spectra from solid state samples. In addition, since electronic devices often operate at high temperatures, it is important to evaluate the thermal dependence of CPL characteristics.

To this end, JASCO has developed various accessories for the CPL-300 spectrophotometer. In this presentation, we describe CPL spectra of solid samples obtained at high temperatures using the CPL-300 spectrophotometer. We also introduce a technique referred to as magnetic circularly polarized luminescence (MCPL) spectroscopy that uses a small permanent magnet.

Keywords : CPL; Europium complex; Solid-state CPL measurement; Temperature-dependent CPL measurements; Magnetic circularly polarized luminescence

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
Fig.1 (A) Sample compartment of CPL-300 with installed JKBR-469 pellet holder, and front view of JKBR-469. (B) CPL and fluorescence spectra of $\text{Eu}(\text{facam})_3/\text{KBr}$ pellet at rotation angles of 0° (green), 45° (blue), and 90° (red). (C) g_{lum} , CPL and fluorescence spectra of $\text{Eu}(\text{facam})_3/\text{KBr}$ pellet (green, left vertical axis) and in DMSO solution (blue, right vertical axis).

- 1) N. Taniguchi, K. Nakabayashi, T. Harada, N. Tajima, M. Shizuma, M. Fujiki, and Y. Imai, *Chem. Lett.*, **2015**, 44, 598-600.
- 2) T. Kimoto, T. Amako, N. Tajima, R. Kuroda, M. Fujiki, and Y. Imai, *Asian J. Org. Chem.*, **2013**, 2, 404-410.
- 3) Y. Okazaki, T. Goto, R. Sakaguchi, Y. Kuwahara, M. Takafuji, R. Oda, and H. Ihara, *Chem. Lett.*, **2016**, 45, 448-450.