

Analysis of membrane protein-specific lipids using gold nanoparticle-based method


(¹Graduate School of Science, Kyushu University) ○Supakorn Wangamnuayporn,¹ Masanao Kinoshita¹, Takayuki Kawai¹, Nobuaki Matsumori¹

Keywords: Gold Nanoparticle, Membrane Protein, Bacteriorhodopsin, Lipid, Liquid chromatography–mass spectrometry

The interaction between membrane proteins (MPs) and surrounding lipids is important to regulate biological properties and stabilize the structure of MPs. In order to understand the MPs-lipid interaction, the analytical method was recently developed based on the surface plasmon resonance (SPR).¹ Briefly, the gold surface of SPR sensor chips was coated by self-assembled monolayer (SAM) containing mercaptohexanol and mercaptohexanoic acid which significantly enhance MP immobilization due to partly hydrophobic environment. After that, MP-lipid interaction was detected by introducing lipid into the MP-immobilized sensor chips. Applying this method to seven-transmembrane bacteriorhodopsin (bR), we found that a halobacterial lipid S-TGA-1 has the highest affinity to bR.² However, lipid isolation from the lipid mixtures extracted from biological membranes is necessary prior to the SPR analysis. Therefore, to overcome the limitation, a new method to identify MP-specific lipids without lipid purification was developed by using SAM-coated gold nanoparticles (AuNPs) based on the idea of SAM-based SPR method.

In this study, AuNPs were coated by SAM, followed by immobilization of bR through the EDC/NHS activation. The effects of SAM chain length, AuNPs size, and the ratio of mercaptoalcohol/mercaptocarboxylic acid were examined for the highest bR immobilization. Then, the lipid mixture extracted from purple membrane was incubated with the bR-immobilized AuNPs. The MP-bound lipid species were quantified by LC-MS and the affinity of each lipid to bR was roughly estimated by introducing the binding index which is the comparison of the MS peak intensity ratio between MP-bound lipids and the extracted lipid mixture. The bR-immobilized AuNPs successfully detected S-TGA-1 as bR-specific lipid which confirmed the result from the SPR-based method. More interestingly, the AuNPs-based method also identified S-TeGA-1 as the lipid with high affinity to bR, providing a new information about the bR-lipid interaction.

In conclusion, the AuNPs-based method, which enables simultaneous lipid screening, does not only reproduce the result for the bR-specific lipid, but also identify a new specific lipid otherwise undetected. This method can be further applied to various MPs to improve understanding on MP-lipid interactions, and thus provide a breakthrough for the study of lipids.



- 1) M. Inada, M. Kinoshita, A. Sumino, S. Oiki, N. Matsumori, *Anal. Chim. Acta* **2019**, 1059, 103-112.
- 2) M. Inada, M. Kinoshita, N. Matsumori, *ACS Chem Biol.* **2020**, 15(1), 197-204.